

Please enter these **calibration parameters** and the **Lot No.** into the BioLecture software!

**pH calibration parameters Lot No. 1824 (BioLector® II/Pro, filter module ID-202)**

Temperature	20°C	21°C	22°C	23°C	24°C	25°C	26°C
φ min	64.51	64.45	64.39	64.34	64.28	64.22	64.16
φ max	17.08	17.08	17.08	17.07	17.07	17.07	17.07
dpH	0.51	0.51	0.51	0.51	0.51	0.51	0.51
pH <sub>0</sub>	6.55	6.54	6.54	6.53	6.52	6.51	6.50
Temperature	27°C	28°C	29°C	30°C	31°C	32°C	33°C
φ min	64.10	64.05	63.99	63.93	63.87	63.81	63.75
φ max	17.07	17.06	17.06	17.06	17.06	17.05	17.05
dpH	0.51	0.51	0.51	0.51	0.51	0.51	0.51
pH <sub>0</sub>	6.50	6.49	6.48	6.47	6.46	6.45	6.45
Temperature	34°C	35°C	36°C	37°C	38°C	39°C	40°C
φ min	63.70	63.64	63.58	63.52	63.46	63.40	63.35
φ max	17.05	17.05	17.04	17.04	17.04	17.04	17.03
dpH	0.51	0.51	0.51	0.51	0.51	0.51	0.51
pH <sub>0</sub>	6.44	6.43	6.42	6.41	6.41	6.40	6.39

**pH sensor properties**

Dynamic range	pH 3.90 - 8.55
Resolution	Up to 0.01 pH (software)
Accuracy	± 0.25 pH at pH 4.35 - 4.95; ± 0.1 pH at pH 4.95 – 7.55; ± 0.25 pH at pH 7.55 - 8.10 (batch calibration)
Response time (t90)	At 25 °C < 30 s
Drift at pH = 7	< 0.005 pH per day (sampling interval of 6 min)
Temperature range	5 °C to 50 °C
Compatibility	Aqueous solutions, ethanol, methanol (max. 5 % v/v)
Sensor stability	sensor material can be degraded by some microorganisms
Cross-sensitivity	Reduced to ionic strength (salinity); high concentration of fluorescent molecules in the visible range can interfere (GFP, (e)YFP); complex media can cause a pH-shift (peptone, yeast extract)
Basic material	pH sensor HP8-1427-02_4 (at least stable for 7 days with CertiPUR-buffer) <b>pH sensors are light-sensitive; please protect them from direct light!</b>

**pH calibration**

Buffer	CertiPUR Reference Material Buffer solutions Set (pH 3.00 ± 0.01 / pH 4.00 ± 0.015 / pH 9.00 ± 0.01 / pH 10.00 ± 0.03, 20 °C); 150 mM Na-Phosphate buffer (16 solutions)
Settings	BioLector protocol = pH-DO-calibration, T = 20-40 °C, 800 rpm, 1000 µL/well, shaking diameter 3 mm, MTP-type = FlowerPlate (MTP-48-BOH)
Calibration device	Hardware ID: BL-02-000F-0032
Calibration phase offset	pH -1.40 (pH Ser.3111-hc, gain 7)
Date of calibration	2018/07/17

**EUROPE**

m2p-labs GmbH  
 Arnold-Sommerfeld-Ring 2 | 52499 Baesweiler | Germany  
 Phone +49-2401-805-330 | Fax: +49-2401-805-333  
 info@m2p-labs.com | support@m2p-labs.com

**USA**

m2p-labs, Inc.  
 400 Oser Ave, Suite 1650 | Hauppauge, NY 11788 | USA  
 Phone +1-631-501-1878 | Fax +1-631-501-1060  
 infoUS@m2p-labs.com | supportUS@m2p-labs.com

Please enter these **calibration parameters** and the **Lot No.** into the BioLecture software!

### DO calibration parameters Lot No. 1824 (BioLector® II/Pro, filter module ID-203)

Temperature	20°C	21°C	22°C	23°C	24°C	25°C	26°C
ϕ cal0	71.84	71.81	71.78	71.75	71.72	71.69	71.66
ϕ cal100	44.69	44.49	44.29	44.09	43.90	43.70	43.50
Temperature	27°C	28°C	29°C	30°C	31°C	32°C	33°C
ϕ cal0	71.63	71.60	71.57	71.54	71.51	71.48	71.45
ϕ cal100	43.31	43.11	42.91	42.71	42.52	42.32	42.12
Temperature	34°C	35°C	36°C	37°C	38°C	39°C	40°C
ϕ cal0	71.42	71.39	71.36	71.33	71.30	71.27	71.24
ϕ cal100	41.93	41.73	41.53	41.33	41.14	40.94	40.74

### DO sensor properties

Dynamic range	0 - 100 % air saturation (a.s.)
Resolution	Up to 0.1 % O <sub>2</sub> (software)
Accuracy	± 5% dissolved oxygen (batch calibration)
Drift at 0% oxygen	< 0.5% O <sub>2</sub> per day (sampling interval of 6 min)
Response time (t90)	< 30 s
Temperature range	5 – 50°C
Sensor stability	sensor material can be degraded by some microorganisms
Cross-sensitivity to	Organic solvents, such as acetone, toluene, chloroform or methylene chloride, Chlorine gas; high concentration of fluorescent molecules in the visible range can interfere (mCherry, tdTomato, dsRed, Nile red); complex media can cause a DO-shift
Basic material	Oxygen sensor Pst3-HG-1742-01 (at least stable for 7 days with CertiPUR-buffer) <b>DO sensors are light-sensitive; please protect them from direct light!</b>

### DO calibration

Calibration	0.5 M Sulfite system (Two-point calibration with oxygen-free environment (sodium sulfite) and air-saturated environment)
Settings	BioLector protocol = pH-DO-calibration, T = 20-40 °C, 800 rpm, 1000 µL/well, shaking diameter 3 mm, MTP-type = FlowerPlate (MTP-48-BOH)
Calibration device	Hardware ID: BL-02-000F-0032
Calibration phase offset	DO -360.25 (DO Ser.4103-hc, gain 7)
Date of calibration	2018/07/17

### Sterilization procedure

Sterilization	Beta irradiation (20 kGy)
BGS-certificate No	518342
Date of sterilization	2018/07/01

#### EUROPE

m2p-labs GmbH  
 Arnold-Sommerfeld-Ring 2 | 52499 Baesweiler | Germany  
 Phone +49-2401-805-330 | Fax: +49-2401-805-333  
 info@m2p-labs.com | support@m2p-labs.com

#### USA

m2p-labs, Inc.  
 400 Oser Ave, Suite 1650 | Hauppauge, NY 11788 | USA  
 Phone +1-631-501-1878 | Fax +1-631-501-1060  
 infoUS@m2p-labs.com | supportUS@m2p-labs.com