

Please enter these **calibration parameters** and the **Lot No.** into the BioLector software!

pH calibration parameters Lot No.2109121 (BioLector® Pro, filter module ID-202/-402)

Temperature	20°C	21°C	22°C	23°C	24°C	25°C	26°C
φ min	62.15	62.11	62.06	62.01	61.97	61.92	61.88
φ max	12.86	12.86	12.85	12.84	12.84	12.83	12.82
dpH	0.58	0.58	0.58	0.58	0.58	0.58	0.58
pH ₀	6.36	6.35	6.35	6.34	6.33	6.32	6.31

Temperature	27°C	28°C	29°C	30°C	31°C	32°C	33°C
φ min	61.83	61.79	61.74	61.69	61.65	61.60	61.56
φ max	12.82	12.81	12.80	12.80	12.79	12.79	12.78
dpH	0.58	0.58	0.58	0.58	0.58	0.58	0.58
pH ₀	6.30	6.29	6.29	6.28	6.27	6.26	6.25

Temperature	34°C	35°C	36°C	37°C	38°C	39°C	40°C
φ min	61.51	61.46	61.42	61.37	61.33	61.28	61.23
φ max	12.77	12.77	12.76	12.75	12.75	12.74	12.73
dpH	0.58	0.58	0.58	0.58	0.58	0.58	0.58
pH ₀	6.24	6.23	6.23	6.22	6.21	6.20	6.19

pH sensor properties

Dynamic range	pH 4.30 - 8.00
Resolution	Up to 0.01 pH (software)
Accuracy	± 0.25 pH at pH 4.80-5.15; ± 0.1 pH at pH 5.15-7.15; ± 0.25 pH at pH 7.15-7.50 batch calibration
Response time (t90)	At 25 °C < 30 s
Drift at pH = 7	< 0.005 pH per day (sampling interval of 6 min)
Temperature range	5 °C to 50 °C
Compatibility	Aqueous solutions, ethanol, methanol (max. 5 % v/v)
Sensor stability	Sensor material can be degraded by some microorganisms
Cross-sensitivity	Reduced to ionic strength (salinity); high concentration of fluorescent molecules in the visible range can interfere (GFP, (e)YFP); complex media can cause a pH-shift (peptone, yeast extract)
Basic material	pH sensor HP8-1811-01_4 (at least stable for 7 days with CertiPUR-buffer) pH sensors are light-sensitive; please protect them from direct light!

pH calibration

Buffer	CertiPUR Reference Material Buffer solutions Set (pH 3.00 ± 0.01 / pH 4.00 ± 0.015 / pH 9.00 ± 0.01 / pH 10.00 ± 0.03, 20 °C); 150 mM Na-Phosphate buffer (16 solutions)
Settings	BioLector protocol = HP8-PSt3-calibration, T = 20-40 °C, 800 rpm, 1000 µL/well, shaking diameter 3 mm, MTP-type = Microfluidic Flower Plate (MTP-MF32C-BOH1)
Calibration device	Hardware ID: BL-09-000F-0032
Calibration phase offset	pH -1.46 (pH Ser. 3111, gain 7)
Date of calibration	2021-08-04

HEADQUARTERS EUROPE

m2p-labs GmbH
Arnold-Sommerfeld-Ring 2
52499 Baesweiler, Germany
Phone +49 -2401 805 330
Fax +49 -2401 805 33
info@m2p-labs.com

SUPPORT

EUROPE
Phone +49 - 2401 805 335
support@m2p-labs.com

N./S. AMERICAS
Phone +1 631 501 1878
supportUS@m2p-labs.com

ASIA PACIFIC
Phone +852 6092 6778
supportAsia@m2p-labs.com

Please enter these **calibration parameters** and the **Lot No.** into the BioLecton software!

DO calibration parameters Lot No.2109121 (BioLector® Pro, filter module ID-203/-403)

Temperature	20°C	21°C	22°C	23°C	24°C	25°C	26°C
ϕ cal0	72.84	72.80	72.76	72.72	72.68	72.65	72.61
ϕ cal100	42.91	42.71	42.50	42.30	42.09	41.89	41.68

Temperature	27°C	28°C	29°C	30°C	31°C	32°C	33°C
ϕ cal0	72.57	72.53	72.49	72.45	72.41	72.37	72.34
ϕ cal100	41.48	41.27	41.07	40.86	40.66	40.45	40.25

Temperature	34°C	35°C	36°C	37°C	38°C	39°C	40°C
ϕ cal0	72.30	72.26	72.22	72.18	72.14	72.10	72.06
ϕ cal100	40.04	39.84	39.63	39.43	39.22	39.02	38.81

DO sensor properties

Dynamic range	0 - 100 % air saturation (a.s.)
Resolution	Up to 0.1 % O ₂ (software)
Accuracy	± 5% dissolved oxygen (batch calibration)
Drift at 0% oxygen	< 0.5% O ₂ per day (sampling interval of 6 min)
Response time (t ₉₀)	< 30 s
Temperature range	5 – 50°C
Sensor stability	Sensor material can be degraded by some microorganisms
Cross-sensitivity to	Organic solvents, such as acetone, toluene, chloroform or methylene chloride, Chlorine gas; high concentration of fluorescent molecules in the visible range can interfere (mCherry, tdTomato, dsRed, Nile red); complex media can cause a DO-shift
Basic material	Oxygen sensor PSt3-HG-1810-01_3 (at least stable for 7 days with CertiPUR-buffer) DO sensors are light-sensitive; please protect them from direct light!

DO calibration

Calibration	0.5 M Sulfite system (Two-point calibration with oxygen-free environment (sodium sulfite) and air-saturated environment)
Settings	BioLector protocol = HP8-Pst3-calibration, T = 20-40 °C, 800 rpm, 1000 µL/well, shaking diameter 3 mm, MTP-type = Microfluidic Flower Plate (MTP-MF32C-BOH1)
Calibration device	Hardware ID: BL-09-000F-0032
Calibration phase offset	DO -360.26 (DO Ser. 4103, gain 7)
Date of calibration	2021-08-04

Sterilization procedure

Sterilization	Beta irradiation (20 kGy)
BGS-certificate No	911184
Date of sterilization	2021-06-24

HEADQUARTERS EUROPE

m2p-labs GmbH
Arnold-Sommerfeld-Ring 2
52499 Baesweiler, Germany
Phone +49 -2401 805 330
Fax +49 -2401 805 33
info@m2p-labs.com

SUPPORT

EUROPE
Phone +49 - 2401 805 335
support@m2p-labs.com

N./S. AMERICAS
Phone +1 631 501 1878
supportUS@m2p-labs.com

ASIA PACIFIC
Phone +852 6092 6778
supportAsia@m2p-labs.com